September 1980

# LABORATION **T6**

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.....  
.....  
.....  
.....

Table 1  
TYPICAL EMISSIVITIES

<u>METALS</u>		<u>FINISHES</u>	
Aluminum, oxidized	.60	Aluminum Paint	.55
Brass, polished	.50	Bronze	.80
Brass, oxidized	.60	Gloss Paint, any color	.90
Copper, oxidized	.65	Dull Paint, any color	.95
Iron, Cast, polished	.25		
Iron, Cast, Oxidized	.70	<u>OTHER</u>	
Iron, wrought	.94	Asbestos	.94
Steel, oxidized	.75	Brick	.90
Steel, Plate, rough	.95	Lampblack	.95
Aluminum, polished	.04	Paper, any color	.95
		Plastics, opaque, any color	.90
		Rubber, any color	.95

*Kirchhoff's law. The upper figure absorbs half the incident radiation. (The lower black body absorbs it all.) The upper figure thus emits half what the black body emits*

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## GENERAL DESCRIPTION

## 1-1 THE SYSTEM

The PR200 radiometer system is built around the pyroelectric detector. The pyroelectric effect occurs in certain classes of materials called ferroelectrics, which possess a permanent electric polarization that is highly dependent on temperature. A small variation of temperature in the detector material produces an electrical current proportional to the rate of temperature change. By mechanically chopping the radiation incident upon the detector an AC current is generated with the same frequency as that of the chopper blade and with an amplitude proportional to the incoming radiant flux. The detector output is applied to a transimpedance amplifier which generates a voltage proportional to the detector current. This combination results in a probe responsivity of approximately 100 V/watt. In the control unit the voltage signal is amplified or attenuated as required and is synchronously detected for optimum signal-to-noise ratio.

## 1-2 PR20 DETECTOR PROBE

The PR20 detector probe contains the detector and its transimpedance amplifier, the chopper motor and its associated control circuitry, the synchronous signal generator, and the calibration L.E.D. A slot in the top of the probe housing permits the interchange of slide-in optical window assemblies. An Infrasil Quartz window assembly is provided with each probe for system calibration and for operation over the .23 to 3.2  $\mu\text{m}$  spectral range. A blank window holder is also provided for wideband windowless operation or custom window mounting. Other types of windows as described in Fig. 6-6 are available from the factory.

The front bezel of the probe rotates to three detented positions spaced 90 degrees apart. In the ZERO position a reference plate, in front of the detector, blocks incoming radiation, allowing the system to be zeroed for internal calibration and absolute measurements. Rotation to the CAL position exposes the detector to the calibration LED for system self test and window transmission compensation. The OPEN position is used for normal external measurements.

## Molelectron Model PR200

## 1-3 PR200 CONTROL UNIT

The PR200 control unit provides power to the PR20 probe and performs all signal processing and control functions. The front panel features both a 3½-digit display for accurate, convenient readout of radiant power and an analog meter for easy adjustment of background null, source intensity, and optical alignment. The NULL control is used for internal zero setting as well as external background cancellation for both hot and cold background temperatures. The range switch provides six decades of system sensitivity.

The WATTS / WATTS/CM<sup>2</sup> switch provides for direct readout of either total radiant flux or irradiance. System response speed is selected by the RESPONSE switch. The FAST position individually optimizes the response speed for each range of sensitivity. The SLOW position provides 100-second integration time for noise reduction on all ranges.

The REAR PANEL includes the probe input connector, BCD output connector labeled DPM output, 0-200 mV full scale analog output jack for data recording, high range (20 and 200 mW) meter ZERO control, gain adjustment control for filter compensation, calibration LED ON-OFF switch, FAC-FIL CAL switch, and power cord connection.

## 1-4 OPTIONAL OPTICAL WINDOWS

Fig. 6-6 shows transmission curves of window assemblies available from Molelectron for use with the PR20.

## 1-5 PR200 SYSTEM SPECIFICATIONS

Active Detector Size:	0.2 cm <sup>2</sup> (5.1 mm dia.)
Detector Position:	5.6 cm behind front bezel edge
Spectral Range and Uniformity:	
Quartz Window:	0.23 to 3.2 μm; ±2%*
Windowless:	0.2 to 20 μm; ±3%*
	0.2 to 40 μm; ±5%*
(See Accessory List for other optional windows)	
Surface Uniformity: (Across Active Area)	±5%

## Molelectron Model PR200

Field of View:	0.1 Sr	
Front Panel Display: Digital Panel Meter:	3½-digit Bipolar LED display providing direct reading of radiant flux in watts and irradiance in watts/cm <sup>2</sup> .	
Analog Indicator Meter:	Provides convenient relative indication for ease of source adjustment, optical alignment, and background nulling.	
Range Selection:	<p>Radiant Flux (Front Panel Switch in WATTS position): 2 μW to 200 mW full scale. Radiation must fall <u>within</u> detector area without exceeding 50 W/cm<sup>2</sup> damage threshold.</p> <p>Irradiance (Front Panel Switch in WATTS/CM<sup>2</sup> position): 2 μW/cm<sup>2</sup> to 200 mW/cm<sup>2</sup> full scale. Radiation must <u>fill</u> entire detector area.</p> <p>Radiance/"Brightness" (Front Panel Switch in WATTS/CM<sup>2</sup> position): 20 μW/cm<sup>2</sup> Sr. to 2 W/cm<sup>2</sup> Sr. full scale. Obtained by multiplying meter reading by factor of 10 with radiation filling entire detector and source filling entire 0.1 steradian field of view.</p>	
Null Adjustment:	+25 μW/cm <sup>2</sup> range for offset of hot or cold backgrounds.	
System Response Time:	Front panel switch selects fast or slow response times (0-90% of full scale) as follows:	
	<u>Range</u>	<u>Fast Response</u> <u>Slow Response</u>
	2 μW	10 sec                  100 sec
	20 μW	3 sec                    100 sec
	200 μW	2 sec                    100 sec
	2,20,200 mW	1 sec                    100 sec
Phase Control:	Screwdriver electronic adjustment recessed on probe rear panel to correct signal for off-axis radiation incident on probe aperture.	

## Molelectron Model PR200

Quartz Window - Typical Lab Background:	Better than $5 \times 10^{-7}$ W/cm <sup>2</sup> per hour after 15 min. warmup.
Electronics Only:	Better than $2 \times 10^{-8}$ W/cm <sup>2</sup> per hour after 15 min. warmup.
System Calibration Accuracy: Windowless or Quartz:	+5% of full scale (with Accessory Windows calibration accuracy depends on spectral uniformity of window transmission)
Linearity:	within 2% up to 50 mW/cm <sup>2</sup> within 5% up to 200 mW/cm <sup>2</sup>
Overall Temperature Coefficient:	Better than 0.2%/°C
Absolute Maximum Radiant Flux (Overload Condition):	10 W
Absolute Maximum Irradiance (Any portion of Active Surface):	50 W/cm <sup>2</sup>
Input Power Requirement:	25 Watts at 115 or 230 Vac, single phase, 50/60 Hz, $\pm$ 10%
Probe Dimensions:	15.0 cm long x 8.2 cm dia (excluding stand)
Control Unit Dimensions:	20.4 cm deep x 17.8 cm wide x 7.6 cm high (excluding stand)
Weight:	Control Unit with Stand, 5 lbs. Probe with Stand, 4 lbs.

## Notes:

\*Spectral Uniformity estimated theoretically based on known characteristics of absorbing coating.

Irradiance measurements valid only if entire active detector surface is irradiated.

Chart recordings of the PR200 output in ZERO position, with no window and with .23 to 35  $\mu$ m optical bandwidth are shown in Figures 6-5A and 6-5C which clearly illustrate the effect of variations in long wavelength energy on the output of the radiometer.

The selection of an optical window that filters out energy at wavelengths that are not required will substantially reduce the effects of thermal drift as indicated by Fig. 6-5B, Strip Chart recording of a PR200 with a quartz window in open air.

## Molelectron Model PR200

## SECTION III

## OPERATION

## 3-1 INTRODUCTION

The wavelength range of the instrument is dependent on the material used for the window. Quartz is supplied with the standard probe for accurate measurements from .23 to 3.2  $\mu\text{m}$ . Measurements outside of the above range may be made with an optional slide-in window assembly or using the detector probe with the blank holder provided for wideband windowless measurements.

## 3-2 INITIAL TURN-ON

Connect the probe to the control unit, plug in the power cord and turn the ON/OFF switch to ON. The digital display should become illuminated. If the display is not illuminated, see Section V for instructions on troubleshooting.

## 3-4 INTERNAL ZERO FOR CALIBRATION AND ABSOLUTE MEASUREMENTS

The NULL control on the control unit front panel adjusts the amplitude and polarity of an electronic signal that is similar to that produced by the detector. Adjustment of the NULL control will either add this signal to or subtract it from the detector signal. In this manner it is possible to NULL out hot or cold background irradiance reaching the detector up to a level of 25  $\mu\text{W}/\text{cm}^2$ .

To zero the instrument internally the front bezel of the probe is rotated to the zero position, as shown in Fig. 6-2. This places a reference plate in the field of view of the detector. The reference plate is the same color temperature as the chopper blade which results in a signal of minimum amplitude. It is this signal that determines the system zero. With the front bezel in the ZERO position adjust the NULL control for a zero reading on the front panel meters. It is helpful to start on the 200  $\mu\text{W}$  range and work down to the 2  $\mu\text{W}$  range. This process is facilitated by observing first the analog meter to determine a coarse zero and then

## Molelectron Model PR200

the digital meter for the fine zero adjustment. In many cases a coarse zero alone will be adequate. The operator should also switch to the FAST mode for the zeroing operation. System zero will not be changed when switched back to SLOW.

If a slide-in window is changed, the system must be re-zeroed before calibration or absolute measurements.

Meter zero for the two high ranges (20 mW and 200 mW) is accomplished by normally zeroing with the NULL control on the 2  $\mu$ W range, then switching to the 20 mW range and zeroing with the rear panel ZERO potentiometer.

## 3-5 INTERNAL SELF TEST

Each PR200 instrument is factory calibrated with NBS traceability. This factory calibration is maintained when the FAC-FIL CAL rear-panel switch is in the FAC (up) position. In this position, the rear panel CAL control potentiometer has no effect. Recalibration by Molelectron or by a competent standards laboratory is recommended every six months, or prior to any particularly critical measurement. A procedure for this calibration is available from Molelectron.

An internal Light Emitting Diode (LED) calibrator provides a useful complete system self-check in the field. The LED calibrator also aids in determining compensation for window or filter transmission.

With the FAC-FIL CAL rear-panel switch in the FAC (up) position, perform the following steps:

- (1) Insert the quartz optical window provided with the radiometer when calibrating the unit.
- (2) After allowing the radiometer to warm up for a minimum of 2 hours with the quartz window in place, zero the instrument as described in Section 3-4.
- (3) Rotate the front bezel of the probe from ZERO to the CAL position.
- (4) Turn on the CAL Light Emitting Diode by switching the CAL switch on the rear panel of the control unit to the ON (up) position. See Fig. 6-3.
- (5) Refer to the calibration label on the probe case. On this label is listed the calibration factor for each window as well as for windowless operation. Note the calibration factor for the quartz window and select the front panel switch positions for the range and units (WATTS or WATTS/CM<sup>2</sup>) as indicated on the label.

## Molelectron Model PR200

- (6) If the display reading agrees with the probe label calibration factor, the instrument is functioning properly. However, a discrepancy of +3% or more indicates need for recalibration by Molelectron or by a competent standards laboratory. Gross discrepancy indicates an instrument malfunction requiring service.
- (7) Turn off the CAL LED by switching the lower CAL switch to the OFF (down) position.

## 3-7-1 RADIANT FLUX (POWER) MEASUREMENTS

Allow the instrument to warm up, calibrate with the desired window (see previous section) switch to the WATTS position, and null the display on the desired range in the absence of the source to be measured.

Make certain that all the radiant flux to be measured falls within the 5 mm diameter detector area and that the source is centered on the probe optical axis and lies within the 0.1 Steradian field of view. The graph in Fig. 6-7 shows the variation in output signal as a point source is moved across the field of view, in both the vertical and horizontal planes.

Irradiate the probe detector with the source radiation and read the radiant flux on the display. Alignment and centering can be checked by gently moving the probe to maximize the display reading. A final correction for off-axis radiation is made by adjusting the PHASE CONTROL\* by means of a small screwdriver through the hole in the probe rear plate until maximum display reading is reached. On the lowest ranges, noise fluctuations may be reduced by switching to the SLOW response speed.

For improved confidence, block or remove the source, readjust the NULL control, and repeat the measurement process. An average of several such readings is recommended.

Do not exceed the maximum display reading of 199.9 mW. Do not exceed a local irradiance of 50 W/cm<sup>2</sup> over any portion of the detector surface.

## 3-7-2 IRRADIANCE MEASUREMENTS

Follow the same procedure as in Radiant Flux Measurements, except for switching to the WATTS/CM<sup>2</sup> position.

In order for irradiance measurements to be valid, it is important that the entire surface of the probe detector be uniformly irradiated. Otherwise the instrument reading will be the average irradiance over the surface area (simply the total radiant flux falling on the detector divided by 0.2 cm<sup>2</sup> area). In general, irradiance measurements are best made with a radiant flux distribution that fills an area larger than the detector, whereas radiant flux measurements are best made where the distribution is entirely within the detector diameter.

For uniform irradiance, do not exceed the maximum display reading of 199.9 mV/cm<sup>2</sup>. For non-uniform irradiance, do not exceed a local irradiance of 50 W/cm<sup>2</sup> at any portion of the detector surface.

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